


CHIROBLOC®

3 models with different ergonomies

3 lead hands without lead.

Made out of a composite bio-compatible material.



CHIROBLOC® 3 - ref. CB3


Weight: 0.27 kg
Warranty: 1 year

Supplied with 5 ForgetMeNot®
see accessories




CHIROBLOC® Large - ref. CBM

Weight: 0.5 kg
Warranty: 1 year



CHIROBLOC® Small - ref. CBP

Weight: 0.3 kg
Warranty: 1 an



Design: e-monceaux.fr / 254 éd.01. © AREX® - 2017

DISTRIBUTED BY




Manufacturer:
AREX®
BP20
3 allée du Clos Tonnerre
F91125 Palaiseau cedex FRANCE

www.arex.fr

www.arex.fr

CHIROBLOC® 3

Made out of a composite bio-compatible material.
Radio-lucency without lead weight.
Will not stain surgical gloves.



Wrist strap /
child positioning

Wrist strap /
adult positioning

CHIROBLOC® Ballast (ref. CBL) : see accessories



One wrist strap: two possible adjustments

5 ForgetMeNot® (FMN2B.Y/ FMN2B.B) 2 functions:

- Reusable finger tourniquet
- ForgetMeNot® units may be used as an attachment to any CHIROBLOC®

ACCESSORIES



CHIROBLOC® Ballast - ref. CBL

It is an independent spare part made out of stainless steel. It may be attached to CHIROBLOC® to weigh it.
Weight: 0.725 kg
Warranty: 5 years



ForgetMeNot®

Yellow color - ref. FMN2B.Y
Blue color - ref. FMN2B.B
Sold by box of 10 units.

FMN is a reusable finger tourniquet, it is sold non sterile.
FMN may also be used as an attachment to any CHIROBLOC®
Warranty: 1 year



ForgetMeNot® is specially useful to replace elastic bands made with surgical gloves.