



Body Aesthetics

Gluteus Implants

Testicular Implants



Calf Implants

Nasal & Malar Implants

Chin Implants

Gluteus Implants

- ✓ 4 different profiles with 34 different variations.
- ✓ Round or anatomical shape.
- ✓ High Cohesive Gel.
- ✓ Textured or smooth shell surface.



ES 62 Smooth Surface

Product Code	Diameter (cm)	Projection (cm)	Volume (cc)
62140	8.4	4.1	140
62180	9.1	4.3	180
62200	9.4	4.4	200
62220	9.7	4.5	220
62240	10.3	4.7	240
62260	10.4	4.7	260
62280	10.6	4.9	280
62300	10.8	50	300

ES 60 Smooth Surface

Product Code	Height (cm)	Width (cm)	Proj. (cm)	Volume (cc)
60150	13.1	9.6	2.3	150
60210	14.6	10.7	2.6	210
60230	12.0	11.2	3.5	230
60275	12.4	11.4	3.8	275
60300	13.9	10.9	3.6	300
60325	14.7	11.3	3.7	325
60375	14.7	11.7	4.1	375
60425	15.0	12.4	4.5	425
60480	15.5	13.2	4.8	480

ES 62C Textured Surface


Product Code	Diameter (cm)	Projection (cm)	Volume (cc)
62140C	8.4	4.1	140
62180C	9.1	4.3	180
62200C	9.4	4.4	200
62220C	9.7	4.5	220
62240C	10.3	4.7	240
62260C	10.4	4.7	260
62280C	10.6	4.9	280
62300C	10.8	50	300

ES 60C Textured Surface

Product Code	Height (cm)	Width (cm)	Proj. (cm)	Volume (cc)
60150C	13.1	9.6	2.3	150
60210C	14.6	10.7	2.6	210
60230C	12.0	11.2	3.5	230
60275C	12.4	11.4	3.8	275
60300C	13.9	10.9	3.6	300
60325C	14.7	11.3	3.7	325
60375C	14.7	11.7	4.1	375
60425C	15.0	12.4	4.5	425
60480C	15.5	13.2	4.8	480

Testicular Implants

- ✓ 5 different variations
- ✓ Oval shape
- ✓ High Cohesive Gel
- ✓ Smooth shell surface
- ✓ Fixation tab eases suturing to surrounding tissues



ES 520			
Smooth Surface			
Product Code	Diameter (cm)	Length (cm)	Volume (cc)
52006	2.2	2.8	6
52012	2.8	3.3	12
52018	3.0	4.2	18
52022	3.3	4.7	22
52026	3.3	5.2	26

Calf Implants

- ✓ 2 profiles with 11 variations
- ✓ Symmetric or asymmetric shape
- ✓ High Cohesive Gel
- ✓ Smooth shell surface
- ✓ Dacron disc on the back facilitates implant insertion




ES 54				
Symmetric Shape				
Product Code	Length (cm)	Width (cm)	Proj. (cm)	Volume (cc)
54060	16.8	3.3	2.0	60
54090	17.0	3.7	2.8	90
54110	17.7	5.8	2.1	110
54120	18.9	5.2	2.6	120
54160	19.5	6.9	2.3	160
54180	19.0	5.9	3.1	180
54220	18.7	7.1	3.2	220




ES 54				
Asymmetric Shape				
Product Code	Length (cm)	Width (cm)	Proj. (cm)	Volume (cc)
54135	22.0	4.9	2.3	135
54140	18.6	5.3	2.6	140
54150	13.9	6.4	2.9	150
54165	23.9	5.6	2.3	165

Nasal & Malar Implants

- ✓ 2 profiles with 4 variations each
- ✓ Nasal or Malar shape
- ✓ Medical-Quality silicone elastomer
- ✓ Fixation holes to ease suturing to surrounding tissue
- ✓ Smooth surface




ES 55		Nasal Implants	
Product Code	Length (cm)	Width (cm)	Base (cm)
55001	5.1	0.8	2.5
55002	5.3	0.8	2.8
55003	5.4	0.8	3.3
55004	5.9	0.9	3.5




ES 33		Malar Implants	
Product Code	Length (cm)	Width (cm)	Thickness (cm)
33001	2.4	1.7	0.3
33002	3.0	2.1	0.4
33003	3.4	2.3	0.4
33004	3.8	2.4	0.5

Chin Implants

- ✓ 2 profiles with 7 variations
- ✓ Medical-Quality silicone elastomer
- ✓ Smooth surface



ES 32 / 34 / 36 / 38		Smooth Surface	
Product Code	Length (cm)	Width (cm)	Thickness (cm)
32000	3.4	1.3	0.9
34000	3.9	1.3	0.9
36000	4.2	1.5	0.9
38000	4.5	1.3	0.7



ES 325		Smooth Surface		
Product Code	Length (cm)	Depth (cm)	Height (cm)	Thickness (cm)
32501	6.2	3.0	0.6	0.6
32502	6.5	3.4	0.8	0.8
32503	7.0	3.8	0.9	0.8